



正 本
ORIGINAL


日期 Date: 2017.02.23 收件日期 Date of Receipt: 2017.02.15
 報告編號 Report No.: TFF6B175 數量 Quantity: 1件 報告頁次/頁數 Page Order/Pages: (P1/1) 來文字號 Ref. No.: 空 白
 報告抬頭 Report Title: 大謙科技材料股份有限公司(U0868) 試件類別 Item: 半成品
 地址 Address: 428 台中市西屯區台灣大道四段1270號18樓之15

試驗報告 土城場區
TEST REPORT TUCHENG

試驗項目	試驗結果	試驗方法
合成血液穿透性 壓力: 120 mmHg	1	無
	2	無
	3	無
	4	無
	5	無
	6	無
	7	無
	8	無
	9	無
	10	無
ASTM F1862-2013		

註: 依委託者所提供來樣資料為: dc 奈米微孔薄膜口罩 材料組合

註: 試驗報告僅就委託者之委託事項提供試驗結果, 不對產品合法性做判斷。



附記: 1. 本報告僅對樣品負責, 樣品保留期限為一個月。

Note: This report is only responsible for the submitted sample(s), which will be kept for one month period.

2. 本報告非經本檢測及驗證部書面同意, 不得摘錄複製, 但全部複製除外。

This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.

3. 試驗報告所載事項不得作為公開廣告及商業推銷之用。

The test report should not be used for public advertisement and commercial promotion.



財團法人紡織產業綜合研究所
授權核發人:
Authorized by president of
Taiwan Textile Research Institute

檢測及驗證部 主任 周國村

紡織所檢測及驗證部 新北市土城區承天路六號
Department of Testing and Certification, Taiwan Textile Research Institute
No.6, Chengtian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R. O. C.)
Tel : +886-2-22670321 ext. 7107, 7110
+886-2-22698229 , +886-2-22698230
Fax : +886-2-22675108 , +886-2-22689839



正 本
ORIGINAL

日期 Date: Feb. 23, 2017 收件日期 Date of Receipt: Feb. 15, 2017

試驗報告 TEST REPORT 土城場區 TUCHENG

報告編號 Report No.: TFF6B175 數量 Quantity: 1PC 報告頁次/頁數 Page Order/Pages: (P1/1) 來文字號 Ref. No.: NIL

報告抬頭 Report Title: DACIAN TECHNOLOGY MATERIAL CO., LTD. (U0868) 試件類別 Item: SEMIFINISHED PRODUCT

地址 Address: 8F-15, NO. 1270, SEC. 4, TAIWAN BLVD., XITN DIST., TAICHUNG CITY, TAIWAN (R.O.C.)

Test Items		Test Results	Test Methods
Synthetic Blood Penetration Pressure: 120 mmHg	1	none	ASTM F1862-2013
	2	none	
	3	none	
	4	none	
	5	none	
	6	none	
	7	none	
	8	none	
	9	none	
	10	none	

Note: Sample description is given by the client: The material structure of dc Face Mask

Note: The test report is the test results issued by the testing institution as requested by the consignor, it shall not determine the legitimacy of the product.

附記: 1. 本報告僅對樣品負責, 樣品保留期限為一個月。

Note: This report is only responsible for the submitted sample(s), which will be kept for one month period.

2. 本報告非經本檢測及驗證部書面同意, 不得摘錄複製, 但全部複製除外。

This report cannot be reproduced in any way, except in full context, without the prior approval in writing of this Department of Testing and Certification.

3. 試驗報告所載事項不得作為公開廣告及商業推銷之用。

The test report should not be used for public advertisement and commercial promotion.



財團法人紡織產業綜合研究所

所長 授權核發人:

Authorized by president of
Taiwan Textile Research Institute

Guo-Tsun Jou

Director.
Department of Testing and
Certification

紡織所檢測及驗證部 新北市土城區承天路六號
Department of Testing and Certification, Taiwan Textile Research Institute
No.6, Chengrlian Rd., Tucheng Dist., New Taipei City 23674, Taiwan (R. O. C.)
Tel : +886-2-22670321 ext. 7107, 7110
+886-2-22698229 , +886-2-22698230
Fax : +886-2-22675108 , +886-2-22689839